

## Product datasheet for **MK203899**

### Uqcrq Mouse qPCR Template Standard (NM\_025352)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** NM\_025352

**Synonyms:** 9.5kDa; 1100001F06Rik; 1500040F11Rik; 5830407L17Rik; AA959903; c1502; QP-C; Qpc; Uqcrb

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP217848) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10<sup>7</sup> copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)