

## Product datasheet for **MK203998**

### Syf2 Mouse qPCR Template Standard (NM\_026780)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_026780
Synonyms:	1110018L13Rik; Cbpin; D4Bwg1551e; Gcipip; mp29; Ntc31; p29
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP216455) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)