

Product datasheet for **MK203999**

Ube2j2 Mouse qPCR Template Standard (NM_001039159)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001039159
Synonyms:	1200007B18Rik; 2400008G19Rik; 5730472G04Rik; AL022923; NCUBE-2; Ubc6; Ubc6p
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP217932) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)