

## Product datasheet for **MK204001**

### Grp2 Mouse qPCR Template Standard (NM\_010815)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** NM\_010815

**Synonyms:** Gads; GRAP-2; GRB2L; GRBLG; GrbX; Grf40; GRID; GrpL; Mona; P38

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP205749) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard:  $50 \times 10^7$  copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)