

Product datasheet for **MK204178**

Ndufab1 Mouse qPCR Template Standard (NM_028177)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_028177

Synonyms: 8kDa; 2210401F17Rik; 2310039H15Rik; 2610003B19Rik; 9130423F15Rik; ACP; CI-SDAP; SDAP

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP208671) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)