

Product datasheet for **MK204189**

Reg3g Mouse qPCR Template Standard (NM_011260)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_011260 |
| Synonyms: | AI449515; REG-3-gamma; reg III-gamma |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (MP212552) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)