

## Product datasheet for **MK204195**

### Casp4 Mouse qPCR Template Standard (NM\_007609)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_007609
Synonyms:	Cas; CASP-4; CASP-11; Casp11; Caspa; Caspl; ich-; ich-3
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP201795) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)