

Product datasheet for **MK204231**

Caly Mouse qPCR Template Standard (NM_026769)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_026769
Synonyms:	0710001P07Rik; 1110004A22Rik; Calcyon; Drdlip
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP201600) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

