

Product datasheet for **MK204250**

4930503B20Rik Mouse qPCR Template Standard (NM_029144)

Product data:

Product Type:	qPCR Template Standards
Accession No:	NM_029144
Synonyms:	1700013A07Rik; 1700095G21Rik; 4930483N21Rik; Gm1652; Samd13
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)