

Product datasheet for **MK204277**

Abat Mouse qPCR Template Standard (NM_172961)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_172961

Synonyms: 9630038C02Rik; AI255750; GABA; Gabaat; Gabat; Gm9851; I54; Laibat; X61497

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP200054) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)