

Product datasheet for **MK204291**

Abhd11 Mouse qPCR Template Standard (NM_145215)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_145215
Synonyms:	1110054D16Rik; A630008N09Rik; Wbscr21
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP200168) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)