

Product datasheet for **MK204409**

Iars Mouse qPCR Template Standard (NM_172015)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_172015
Synonyms:	2510016L12Rik; AI327140; AU044614; E430001P04Rik; ILRS
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP206872) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

