

Product datasheet for **MK204410**

Pink1 Mouse qPCR Template Standard (NM_026880)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_026880
Synonyms:	1190006F07Rik; AU042772; AW557854; BRPK; mFLJ00387
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP212207) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)