

Product datasheet for **MK204481**

Sae1 Mouse qPCR Template Standard (NM_019748)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_019748
Synonyms:	2400010M20Rik; 2610044L12Rik; AL033372; AOS1; AW743391; D7Ertd177e; HSPC140; Sua1; Uble1a
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP216380) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)