

Product datasheet for **MK204498**

Pmpca Mouse qPCR Template Standard (NM_173180)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_173180
Synonyms:	1200002L24Rik; 4933435E07Rik; Alpha-MPP; P-55
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP210674) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)