

## Product datasheet for **MK204563**

### Syp Mouse qPCR Template Standard (NM\_009305)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** NM\_009305

**Synonyms:** A230093K24Rik; AI848995; p38; Syn

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP216473) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard:  $50 \times 10^7$  copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)