

Product datasheet for **MK204608**

Atp6ap1 Mouse qPCR Template Standard (NM_018794)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_018794
Synonyms:	16A; AC45; AI316502; Atp6ip1; Atp6s1; AW108110; C7-1; CF2; mFLJ00383; VATPS1; XAP-3
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP201083) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)