

Product datasheet for **MK204637**

Aldh4a1 Mouse qPCR Template Standard (NM_175438)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_175438
Synonyms:	A930035F14Rik; Ahd-1; Ahd1; Aldh4; Aldh5a1; E330022C09; P5cd; P5cdh; P5cdhl; P5cdhs; Ssdh1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP200517) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)