

## Product datasheet for **MK204734**

### Arhgdia Mouse qPCR Template Standard (NM\_133796)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_133796
Synonyms:	5330430M07Rik; Gdi-1; Rho-GDI; RhoDGI; RhoGDI-1
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP212758) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)