

## Product datasheet for **MK204744**

### **Dnaja2 Mouse qPCR Template Standard (NM\_019794)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**ACCN:** NM\_019794

**Synonyms:** 1500017M13Rik; 2010206B19Rik; DNAJ; Dnaj3; DNJ3; HIRIP4; mDj3; PRO3015

**Components:** 1 vial of lyophilized SybGREEN qPCR primer (MP203839) mix (1 nmol each primer, sufficient for 200 reactions)  
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10<sup>7</sup> copies, double-stranded DNA  
1 vial containing 1 mL 10x PCR template standard dilution buffer

**Quality Control:** The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)