

## Product datasheet for **MK204761**

### Sf3b3 Mouse qPCR Template Standard (NM\_133953)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_133953
Synonyms:	1810061H24Rik; 5730409A01Rik; AA409318; D8Ertd633e; mKIAA0017; RSE1; SAP130
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP216157) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)