

Product datasheet for **MK204794**

Ap1m2 Mouse qPCR Template Standard (NM_009678)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_009678 |
| Synonyms: | D9Erttd818e; mu1B; [m]1B |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (MP220985) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)