

Product datasheet for **MK204814**

Scarb1 Mouse qPCR Template Standard (NM_016741)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_016741
Synonyms:	AI120173; CD36; Cd36l1; Chohd1; Cla-1; Cla1; D5Ertd460e; Hdlq1; Hlb398; mSR-BI; SR-B1; SR-BI
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP215005) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)