

Product datasheet for **MK204859**

H3c13 Mouse qPCR Template Standard (NM_178215)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_178215

Synonyms: H3-615; H3-616; H3c2; H3c3; H3c4; H3c6; H3c7; H3c14; H3c15; Hist2h; Hist2h3b

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP206231) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)