

## Product datasheet for **MK204890**

### **Eef1a1 Mouse qPCR Template Standard (NM\_010106)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_010106
<b>Synonyms:</b>	MGC7551; MGC8115; MGC8209; MGC18758; MGC27859; MGC102592; MGC103271; MGC118397
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (MP204369) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)