

Product datasheet for **MK204944**

Trem2 Mouse qPCR Template Standard (NM_031254)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_031254

Synonyms: Trem; TREM-2; Trem2a; Trem2b; Trem2c

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP217538) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50 X 10⁷ copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)