

Product datasheet for **MK204966**

Gnpnat1 Mouse qPCR Template Standard (NM_019425)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_019425

Synonyms: AU017428; AU040593; EMeg32; Gpnat1; Gsnpat

Components: 1 vial of lyophilized SybGREEN qPCR primer (MP205491) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)