

Product datasheet for **MK205019**

Dpp9 Mouse qPCR Template Standard (NM_172624)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_172624
Synonyms:	6430584G11Rik; A330078I11; DPP IX; DPRP2
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP203703) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)