

Product datasheet for **MK205096**

Zbtb20 Mouse qPCR Template Standard (NM_019778)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_019778
Synonyms:	1300017A20Rik; 7330412A13Rik; A930017C21Rik; D16Wsu73e; DPZF; HOF; ODA-8S; Oda8; Zfp288
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP218677) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)