

## Product datasheet for **MK205110**

### **Pou5f1 Mouse qPCR Template Standard (NM\_013633)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_013633
<b>Synonyms:</b>	NF-A3; Oct; Oct-; Oct-3; Oct-3/; Oct-3/4; Oct-4; Oct3; Oct3/; Oct3/4; Oct4; Otf; Otf-; Otf-3; Otf-4; Otf3; Otf3-; Otf3-rs7; Otf3g; Otf4
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (MP211358) primer mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)