

Product datasheet for **MK205113**

Myc Mouse qPCR Template Standard (NM_010849)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_010849
Synonyms:	AU016757; bHLHe3; bHLHe39; Myc2; N; Niard; Nird
Components:	1 vial of lyophilized SybGREEN qPCR primer (MP208493) primer mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)