

Product datasheet for **MK300015**

Let7f-1 Mouse qPCR Template Standard (MI0000562)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0000525

Synonyms: Let7f-1 mmu-let-7f-1; MI0000562; MIMAT0000525; mmu-let-7f

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)