

## Product datasheet for **MK300035**

### Mir1193 Mouse qPCR Template Standard (MI0006298)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0005851
Synonyms:	MI0006298; MIMAT0005851; Mir1193 mmu-mir-1193; mmu-miR-1193
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)