

## Product datasheet for **MK300047**

### Mir127 Mouse qPCR Template Standard (MI0000154)

#### Product data:

|                  |                                                                                                                                                                                                       |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:    | qPCR Template Standards                                                                                                                                                                               |
| ACCN:            | MIMAT0004530                                                                                                                                                                                          |
| Synonyms:        | MI0000154; MIMAT0004530; Mir127 mmu-mir-127; mmu-miR-127                                                                                                                                              |
| Components:      | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment                                                                                                     |
| Storage:         | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.                                                                                                   |



[View online »](#)