

## Product datasheet for **MK300084**

### Mir16-1 Mouse qPCR Template Standard (MI0000565)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0000527
Synonyms:	MI0000565; MIMAT0000527; Mir16-1 mmu-mir-16-1; mmu-miR-16
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)