

## Product datasheet for **MK300097**

### Mir188 Mouse qPCR Template Standard (MI0000230)

#### Product data:

|                  |   |
|------------------|---|
| Product Type:    | qPCR Template Standards   |
| ACCN:            | MIMAT0000217  |
| Synonyms:        | MI0000230; MIMAT0000217; Mir188 mmu-mir-188; mmu-miR-188-5p   |
| Components:      | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment   |
| Storage:         | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)