

Product datasheet for **MK300131**

Mir1937a Mouse qPCR Template Standard (MI0009926)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0009401
Synonyms:	MI0009926; MIMAT0009401; Mir1937a mmu-mir-1937a; mmu-miR-1937a
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)