

Product datasheet for **MK300150**

Mir1950 Mouse qPCR Template Standard (MI0009941)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0009417

Synonyms: MI0009941; MIMAT0009417; Mir1950 mmu-mir-1950; mmu-miR-1950

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)