

Product datasheet for **MK300164**

Mir1964 Mouse qPCR Template Standard (MI0009961)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0009437

Synonyms: MI0009961; MIMAT0009437; Mir1964 mmu-mir-1964; mmu-miR-1964

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)