

Product datasheet for **MK300206**

Mir22 Mouse qPCR Template Standard (MI0000570)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0004629

Synonyms: MI0000570; MIMAT0004629; Mir22 mmu-mir-22; mmu-miR-22

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)