

Product datasheet for **MK300234**

Mir299 Mouse qPCR Template Standard (MI0000399)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0000377

Synonyms: MI0000399; MIMAT0000377; Mir299 mmu-mir-299; mmu-miR-299

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)