

Product datasheet for **MK300236**

Mir29b-1 Mouse qPCR Template Standard (MI0000143)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0004523
Synonyms:	MI0000143; MIMAT0004523; Mir29b-1 mmu-mir-29b-1; mmu-miR-29b
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)