

## Product datasheet for **MK300240**

### Mir301a Mouse qPCR Template Standard (MI0000401)

#### Product data:

|                  |                                                                                                                                                                                                       |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:    | qPCR Template Standards                                                                                                                                                                               |
| ACCN:            | MIMAT0000379                                                                                                                                                                                          |
| Synonyms:        | MI0000401; MIMAT0000379; Mir301a mmu-mir-301a; mmu-miR-301a                                                                                                                                           |
| Components:      | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment                                                                                                     |
| Storage:         | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.                                                                                                   |



[View online »](#)