

Product datasheet for **MK300252**

Mir32 Mouse qPCR Template Standard (MI0000691)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0000654

Synonyms: MI0000691; MIMAT0000654; Mir32 mmu-mir-32; mmu-miR-32

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)