

Product datasheet for **MK300280**

Mir351 Mouse qPCR Template Standard (MI0000643)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0000609

Synonyms: MI0000643; MIMAT0000609; Mir351 mmu-mir-351; mmu-miR-351

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)