

Product datasheet for **MK300283**

Mir363 Mouse qPCR Template Standard (MI0000765)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0000708

Synonyms: MI0000765; MIMAT0000708; Mir363 mmu-mir-363; mmu-miR-363

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)