

Product datasheet for **MK300292**

Mir376c Mouse qPCR Template Standard (MI0003533)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0005295

Synonyms: MI0003533; MIMAT0005295; Mir376c mmu-mir-376c; mmu-miR-376c

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)