

Product datasheet for **MK300304**

Mir412 Mouse qPCR Template Standard (MI0001164)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0001094

Synonyms: MI0001164; MIMAT0001094; Mir412 mmu-mir-412; mmu-miR-412

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)