

Product datasheet for **MK300328**

Mir466b-1 Mouse qPCR Template Standard (MI0005502)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0004875

Synonyms: MI0005502; MIMAT0004875; Mir466b-1 mmu-mir-466b-1; mmu-miR-466b-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)