

## Product datasheet for **MK300349**

### Mir470 Mouse qPCR Template Standard (MI0002405)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** MIMAT0002111

**Synonyms:** MI0002405; MIMAT0002111; Mir470 mmu-mir-470; mmu-miR-470

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)